

# **Advanced Biophotonics Sensors for Environmental and Biomedical Applications**

**Tuan Vo-Dinh and Fei Yan**


*Center for Advanced Biomedical Photonics  
Oak Ridge National Laboratory  
Oak Ridge, TN 37831  
U.S.A.*

*E-mail: [vodinht@ornl.gov](mailto:vodinht@ornl.gov)*

# Outline

- **Introduction to Raman spectroscopy**
- **Surface-enhanced Raman scattering gene probes**
- **Raman integrated tunable sensor**
- **Optical nanobiosensors for single cell analysis**
- **Multi-functional biochips**
- **“Optical biopsy” based on laser-induced fluorescence**
- **Conclusion**

# Raman scattering



- **Based on vibrational transitions**
- **Offers structural information**
- **Narrow bands potentially afford spectral selectivity**
- **Widely considered as complement to IR spectroscopy**
- **Continuous, tunable excitation sources are not necessary**
- **Signals generally very low**

# Surface-enhanced Raman Scattering (SERS)

- Produces Raman scattering enhancement of generally  $10^6$ - $10^8$
- Enhancements as high as  $10^{15}$  observed with specific compounds in single-molecule detection studies
- Results from the adsorption of compounds on submicron-textured metal surfaces
- Has been explained in terms of electromagnetic and molecular enhancement models

$$\mu = E \cdot \alpha$$


$\mu$  = Induced dipole moment

$E$  = Incident electric field

$\alpha$  = Molecular transition polarizability

# SEM Photographs of SERS-active nanoparticles


A) 0.33 %




B) 1 %





C) 3.3 %



D) 10 %




# Calculated Raman enhancement and excitation profile using models for solid metal nanosphere and concentric dielectric/metallic nanospheres




Source: M. Kerker, *Acc. Chem. Res.* 17 (1984) 271

# Schematic diagram of a tapered fiberoptic SERS probe with silver island SERS-active layer




# Microparticle-based SERS substrate





**SEM micrograph of silver-coated polystyrene microspheres**

## Microparticle-based Substrate Parameters:


- **Microparticle material** (e.g. alumina, titanium dioxide, polystyrene, fumed silica)
- **Microparticle size** (e.g. 0.05 - 0.5  $\mu\text{m}$ )
- **Metal** (e.g. silver, gold, copper)
- **Metal thickness** (e.g. 750 - 1500  $\text{\AA}$ )



# Single-molecule detection of cresyl fast violet using silver islands SERS substrate



# Instrumental set-up for spectral recording of individual spots of SERGen probes




An individual spot corresponds to an individual microdot on a hybridization platform


Readily assembles using commercially available or off-the-shelf components (10K)

Bandpass filter isolates laser line prior to sample excitation. The collection module includes a Raman holographic filter which rejects the Rayleigh scattered radiation.


# DNA sequence used with SERGen probe for HIV detection




# SERS gene probe development via PCR and hybridization for selective detection




# Demonstration of the SERGen Technique for the Detection of HIV1 gag gene




# SERS Spectrum of BCB




# Simultaneous detection of HIV and HEPC SERS gene probe hybridizations



# SERS spectra of Raman active dyes



# Multiplexed DNA assays using SERS-labeled PCR amplified products



# SERS Challenges

Experimentally:

- Metallic nanostructures with nano-meter precision (nano-patterning)

Theoretically:

- Better understanding of Raman scattering enhancement mechanism

**RAMITS**

*Raman Integrated Tunable Sensor*

# Raman instrument for field analysis

## Requirements for field analysis:


- Provides rapid, on-site, in situ analysis,
- Improve the quality, reduce the costs, and enhance the safety of characterization and monitoring activities (i.e., hazardous waste sites).

## Analytical advantages of Raman technique:

- Produces a unique, spectral "fingerprint" for each analyte
- Analyzes samples such as solids, liquids, and sludges directly, without the need for isolation or cleanup.

# RAMiTS:

*RAMan Integrated Tunable Sensor*



## System Specifications

### Size:

8.0 x 9.0 x 21.0 in

### Weight:

~ 39 lbs

### Resolution:

~ 7.5  $\text{cm}^{-1}$



### Battery Lifetime:

3 hours


### User Interface:

Touchscreen

# Optical System




# AOTF Operation







# Image of the Optical Housing



# Power Control Circuitry




# Raman Spectrum of Benzene, Toluene and Napthalene Mixture





# RAMiTS for the differentiation of isomers

- Capable of rapid analyses, total acquisition for each of the spectra shown below is approximately 3 seconds or 16 ms per point.
- High resolution Raman spectra allow for the analysis of complex mixtures.



# RAMiTS for the detection of pesticides



- Capable of trace detection of **chemical pesticides** using specially designed surface-enhanced Raman (SERS) snap-top vials.

# RAMiTS

## RAMan Integrated Tunable Sensor

- Capable of performing analyses in seconds.
- Automated chemical identification (spectral library).
- Suitable for remote and non-invasive analyses.
- Capable of operating under typical room light conditions with no significant background signal
- Waterproof, complying with NEMA-4 standards




**FAA and dozens of companies are actively pursuing license agreements for the use of or commercializing this device**

# **Optical Nanobiosensors**


*An Enabling Technology for  
Single Cell Analysis*

# Cell –The Basic Unit of Life



- Two ways to perform Single Cell Analysis:
  - (i) Extraction of its component parts from a whole cell, followed by separation and MS analysis
  - (ii) Examination of a whole cell its elements and their relationship
- **ABST-ORNL: developed OPTICAL NANOSENSORS for single living cell analysis**

# Concept of optical nanosensor for single cell analysis



**Analyte**


**Bioreceptor**

**Molecular  
Recognition**

**Transducer**


**Measurement**

**Data Recording  
and Display**




# Steps in nanosensor fabrication


Fabrication of Nanotips ( $600\ \mu\text{m} \rightarrow 50\ \text{nm}$ )



Silver Coating (~ 100-200 nm)



Silanization of Fiber nanotips (GOPS)  
Glycidoxypropyltrimethoxysilane




Activation of Fiber nanotips (CDI)  
Carbonyldiimidazole




Protein Binding (e.g., mouse IgG for BaP)  
Benzo[a]pyrene

# Protein Binding Chemistry




- Functionalize tips to facilitate protein immobilization

# SEM Image of a Pulled Fiber



# How the Optical Nanosensor Works?

- Total internal reflection (TIR) is used to generate an evanescent field in the near field of the optical nanosensor that is used to excite the analyte of interest bound onto the optical nanosensor
- The evanescent field provides a highly localized excitation area in the near field of the optical nanosensor, therefore good S/N ratio




# Various Assay Formats for Single Cell Analysis Using Optical Nanosensors

## (I) Direct Antibody Binding: Chemical Analysis




## (II) Indirect Antibody Binding: Protein Analysis



## (III) Peptide Substrate: Protein Analysis



# Intra-cell Workstation: Nikon Inverted Fluorescence Microscope




NIKON inverted fluorescence microscope



MMO-203 nanosensor holder 3-axis fine movement micromanipulator

# Optical Nanosensor

## *In vivo* measurement of Benzo[a]Pyrene



- Fluorescence signal produced when Anti-BaP/BaP complex is excited is converted to a measurable electrical signal and recorded





Fig 6


# BaP Measurements of Human Mammary Carcinoma Cells



# Single Cell Analysis: Enzyme Peptide Substrate-based Nanosensor for Monitoring Apoptosis


- Apoptosis or programmed cell death is induced in human mammary carcinoma (MCF-7) cells
- Aminolevulinic Acid (ALA), a photodynamic therapy drug is used to induce apoptosis
- Once apoptosis is induced, the near-field of optical nanosensor is used to probe the cytoplasm of individual human mammary carcinoma cells for caspase-9 activity

# Aminolevulinic Acid-induced Apoptosis




ALA photoactivation process in a MCF-7 cell, showing the activation of caspase-9 and ultimately leading to apoptosis.

# Application: Enzyme Substrate-based Optical Nanosensors for Caspase-9 Activity







- Monitoring enzyme activity requires the use of specific fluorescent substrates such as:

- Leu-Glu-His-Asp-7-amino-4-methyl-coumarin (LEHD-AMC) that is specific of Caspase-9
- LEHD-AMC is not fluorescent but the enzyme activity liberates the fluorescent compound, AMC ( $\lambda_{ex}$  341nm,  $\lambda_{em}$  441nm)




# Intracellular Measurement




**LEFT TO RIGHT, single cell analysis, optical nanosensor outside MCF-7 cell (i) right before insertion into MCF-7 cell (ii) inserted into MCF-7 cell and (iii) extracted from MCF-7 cell**

# Mitosis



After probing the cytoplasm of the cell with a nanoprobe...

# Mitosis




The cell carried out normal cellular activities, specifically mitosis

## *In vivo* measurement of apoptosis using tetrapeptide fluorogenic substrate


- 2 groups:
  - Control: Functionalized optical nanosensors with peptide substrate attached are incubated in cells in which apoptosis has not been induced
  - Experimental: Functionalized optical nanosensors are incubated in cells undergoing apoptosis
- ALA was used to induce apoptosis in MCF-7 cells
- Cells were placed on the stage of the microscope and intracellular measurements were made

# Measurement of Caspase-9 activity

[Detection of 7-amino-4-methoxy coumarin (AMC)]




# *In vivo* measurement of Caspase-9 activity



**For the first time, we identified the onset of the mitochondrial pathway of apoptosis in a single live cell by detecting the enzymatic activity of caspase-9.**

# **The Multifunctional Biochip**

# The Biochip Technology





- 2-D array of independently operating photodiodes
- On-board signal amplification and data treatment
- Complementary metal oxide silicon (CMOS)-based microelectronics integrated onto a single platform
- Coupled to compact sampling system

## ADVANTAGES


- Compact design
- Multiple assays possible on single platform
- Increased throughput
- Microscale sampling capability
- Low power consumption
- Cost effectiveness

# Demonstration of multifunctional capability of the integrated circuit (IC) biochip




Bioreceptors: antibody, DNA, enzymes, tissues, and organelles etc.

# Biochip and Miniature Diode Laser




- Biochip mounted on printed circuit board shown with miniature diode laser used for excitation
- Excitation wavelength of diode laser is 635 nm
- Diode laser power is 5 mW and uses a 5 V, 60 mA power supply

# Schematic Diagram of 8 x 8 Biochip With 4 Different Designs of IC Signal Processors




# Schematic Diagram of Design Q2 of the 8 x 8 Biochip



**4 x 4 Photodiode Array**


**IC Processor**

**Integrating Capacitor**




**Operational Amplifier**

# Modified Ink-Jet Printer Technology for Rapid Formatting of Sample Platforms



- The parallel deposition process for multiple capture probes enables very rapid production of complex sample platforms.
- This figure illustrates Biochip signals observed for a DNA-based assay for the FHIT gene, performed on a sample platform which was preformatted using the modified ink-jet printer. The printer was used to produce 20 nl (0.8-1.0 mm dia.) dots of both FHIT and E4 (control) capture probes.

# Sensitive on-chip detection of *E. coli* using the biochip and biofluidic system



- An on-chip assay for *E. coli* has been performed using an antibody probe and 30-min incubation.
- A proportional response to the number of *E. coli* organisms deposited on a preformatted sampling platform has been observed.
- As few as 21 *E. coli* organisms have been detected through on-chip monitoring of a sandwich assay.

# Portable Biochip Device

- Currently licensed to *HealthSpex, Inc.*, a local company in Knoxville, TN.
- Expected to appear in hospital/doctor's office in a couple of years




# **Laser-Induced Fluorescence (LIF) for Cancer Diagnosis**


*OPTICAL BIOPSY*

# BACKGROUND


- Esophageal cancer is one of the most prevalent types of cancer
- In general, endoscopy is used to detect malignancies in the gastro-intestinal (GI) tract
- Once a suspected tumor is founded, it is removed or biopsies are taken for determination of histopathology
- The laboratory results are generally not available for several days




**Reflux of gastric acid, pepsin, and other components of gastric juice into the esophagus, typically through a hernia, causes normal squamous epithelium (*top*) to be replaced by Barrett's metaplastic epithelium (*middle*). Endoscopic biopsy is necessary for the diagnosis of Barrett's esophagus.**



**Systematic endoscopic protocol,** involving biopsies for pathology and for flow cytometry, permits mapping of the Barrett's esophageal segment in patients undergoing surveillance.




The measurement takes 0.6 s for each tissue site




FROM SILVERSTEIN AND TYTGAT, WITH PERMISSION.

**Barrett's esophageal adenocarcinoma with intramural metastasis is visible endoscopically from a proximal position**


Normalized Fluorescence Curves  
(Relative Intensity)



Differential Normalized Fluorescence (DNF)



# N-FILES UPPER G.I. INTENSITY @ 480nm



**Excellent agreement (98%) with histopathology in classification of normal tissue and malignant tumors of GI cancer in clinical studies involving over 100 patients**

**\*Computer-assisted tomography (CAT) scan reveals lung cancer spreading into areas underlining the esophageal mucosa**

# Conclusion

- The SERS technology has great promise for use in pathogen diagnosis, as well as in environmental monitoring.
- The RAMiTS device integrates the optical module, computer and touch screen interface into a compact, user-friendly, free-standing Raman detection unit.
- By specifically developing optical nanosensors for intracellular measurement applications and using methods analogous to those currently used in cellular and molecular biology, the transfer of technology between scientific disciplines is enhanced.

## Conclusion (cont'd)

- An IC Biochip device has been developed which is very compact and cost effective, with potential for high throughput and multifunctional capability
- Optical biopsy based on LIF will hopefully preclude the need for unnecessary tissue biopsy and thus contribute immensely to early Barrett's esophageal adenocarcinoma detection critical for the timely treatment.

*“After years of intensive effort in basic research, spectroscopy is now poised to emerge as a clinically useful tool (Gastroenterology Editor comments)”*

# Advanced Biomedical Science and Technology Group

T. Vo-Dinh, Ph.D., *Group Leader*<sup>a</sup>

J. B. Cooper, *Secretary*

## Advanced Monitoring Development (AMD)

T. Vo-Dinh, Ph.D.  
G. Griffin, Ph.D.  
B. Hingerty, Ph.D.<sup>3</sup>  
A. L. Wintenberg, Ph.D.<sup>4</sup>  
B. M. Cullum, Ph.D.<sup>1</sup>  
J. Mobley, Ph.D.<sup>1</sup>  
J. Song, Ph.D.<sup>1</sup>  
D. L. Stokes, Ph.D.<sup>1</sup>  
D. N. Stratis, Ph.D.<sup>1</sup>  
P. Kasili,<sup>2</sup>  
M. W. Williams, Ph.D.<sup>6</sup>  
M. J. Sepaniak, Ph.D.<sup>6</sup>  
P. Viallet, Ph.D.<sup>6</sup>

## Biosystems Modeling (BS)

K. F. Eckerman, Ph.D.  
C. E. Easterly, Ph.D.  
R. W. Leggett, Ph.D.  
N. B. Munro, Ph.D.<sup>3</sup>

## Forensics (F)

A. A. Vass, Ph.D.

## Ionizing Radiation Dosimetry & Metrology (IRDM)

J. S. Bogard, Ph.D.

## Photophysics, DNA Physics (PDNA)

C. H. Chen, Ph.D.  
S. L. Allman  
L. J. Sammartano, Ph.D.<sup>6</sup>

## Structural Biology (SB)

G. J. Bunick, Ph.D.  
L. Cacheiro<sup>6</sup>  
B. L. Hanson, Ph.D.<sup>6</sup>  
J. M. Harp, Ph.D.<sup>6</sup>

a Corporate Fellow

1 Postdoctoral Associates

2 Graduate Student

3 Part-Time

4 On loan part-time from ESTD

5 On Loan from I&C (50%)

6 Visiting Scientists/Guest Assignments

## INSTRUMENTATION

- Advanced laser techniques
- Biological mass spectrometry
- Calibrated radiation exposure
- Crystallography
- Microbiology techniques
- Nanosensors, biosensors, biochips
- Scattering techniques (light, X-ray, neutron)
- Tissue culture techniques

## THEORETICAL MODELING

- Computational modeling:
  - biomolecules (DNA, proteins)
  - crystallographic structures
- Mathematical modeling:
  - biokinetics and dosimetry
- Tissue-light interactions modeling

## APPLICATIONS

- Biomedical diagnostics
- DNA/Gene expression decoding
- Health standards
- Low-dose exposure
- Pathogen sensing
- Structural biology of complex systems

# Acknowledgments

- **U. S. Department of Energy**
- **National Institutes of Health**
- **U. S. Department of Defense**
- **U. S. Department of Justice**

Questions ?